



Offer your patients KC's most comfortable mammogram.

Clarity Care, your new partner for advanced breast imaging services.



Finally, a More Comfortable Mammogram

Women have often complained about the compression pain felt during their mammogram. The wait is finally over.



At Clarity Care, patient comfort is our top priority. That's why **our 3-D mammography program features Fujifilm's Patented Comfort Paddle**, designed to enhance the mammogram experience with a softer, more even compression.

Comfort Paddle Key Benefits

Designed with your patients in mind, the Comfort Paddle offers a range of features that make our 3D-mammography experience more comfortable and patient-friendly:



Gentler, More Even Pressure: The innovative paddle design distributes compression more evenly across the breast, helping to reduce discomfort compared to traditional flexible paddles.



Adaptive Fit: A slotted edge allows the paddle to flex naturally, contouring from the center of the breast to the thinner edges for a better fit.



More Gradual Compression: The tilting action of the paddle ensures a smoother, more even compression process, improving both comfort and image quality.

Here's what our patients are saying about their mammogram experience at Clarity Care:

“The easiest and **most comfortable mammogram** that I have ever experienced! Easy to make appointment, friendly staff and nice facility. I have already recommended to a friend of mine.”
- Kathy H., Clarity Care Patient



Comfort without Sacrificing Quality

While we value patient experience at Clarity Care, we also prioritize quality and safety. Our breast center uses Fujifilm's ASPIRE Cristalle mammography system, designed to reduce radiation exposure while improving image quality:



Lower Dose, Safer Exams

The system's advance detector design has a special pixel layout that helps capture x-ray signals more efficiently, allowing us to use up to **50–70% less radiation compared to many traditional systems** and depending on patient breast density and imaging mode.*



Personalized Dose Settings

Our technologists can **adjust dose levels based on each patient's needs**, ensuring a safe, tailored experience for every woman we screen.*



Smarter Imaging with AI

The system's automatic image optimization uses intelligent technology to analyze each patient's breast tissue — even if there are implants or dense tissue — and **automatically adjusts the exposure for the best possible images.***



Enhanced Image Clarity

The AI-powered software **fine-tunes images** after they are taken, making important structures like glands, skin lines, and any potential abnormalities clearer and easier to detect, especially in women with dense breast tissue.*



Built-In Quality Checks for Consistency

Every scan goes through an **automated quality control** process to ensure that images are clear, consistent, and ready for an accurate reading — helping us deliver the highest standard of care to your patients.*

* Based on data from the Fujifilm ASPIRE Cristalle White Paper (DOC-0005190-C). Includes dose reductions up to 70%, improved detector sensitivity, exposure accuracy via iAEC, and enhanced image quality with AI-based processing. Performance may vary based on breast density, positioning, and imaging mode.



Convenient & Accessible



Partner with Clarity Care

At Clarity Care, we are committed to delivering precise diagnostics and effective treatments, ensuring optimal patient outcomes, while supporting you and your patients every step of the way. We invite you to consider Clarity Care as an extension of your practice, providing your patients with advanced imaging and minimally invasive procedures in a comfortable doctor's office setting.

For referrals, please call our office at **913.944.4900**, or fax an order to:
Radiology Fax: 913-701-3714
Specialty Procedure Fax: 913-701-3735



To speak with our physician liaison or request to access to our physician portal, please contact **Melissa Guzman** at **913.944.4901** or **mguzman@claritycarekc.com**.